



Enhance traceability

Comprehensive sample identification demands advanced printing capabilities.

Epredia™ VEGA laser cassette printers offer:

- High-resolution, durable printing
- Up to 100-character 2D barcodes
- Maximum laser print speed of 2.2 seconds per cassette printing
- Intelligent cassette selection from six-hopper or eight-hopper carousel
- Quiet printing, < 65 dB





- High-resolution, 600-dpi barcode to include more information
- High-quality barcode printing may reduce instances of reprints and non-scannable barcodes



- LIS connectivity to streamline workflow
- Easy-to-use connection for smooth installation



- Reduce operational cost with no printhead replacements
- Robust laser technology consistently prints within 2.2 seconds for a faster turnaround time

Ordering Information

Epredia Ref.	
Epredia VEGA Laser Cassette Printer (stand alone unit)	
Epredia Laser Cassette Printer with 6 hoppers	8000038
Epredia Laser Cassette Printer with 6 hoppers with Sorter	8000040
Epredia Laser Cassette Printer with 6 hoppers with DUO	8000041
Epredia Laser Cassette Printer with 8 hoppers	8000062
Epredia Laser Cassette Printer with 8 hoppers with Sorter	8000063
Epredia Laser Cassette Printer with 8 hoppers with DUO	8000064
Accessories	
Slide Delivery system for Vega	9001140
Cassette Sorter for Vega	8000042
Cassette collector for Vega Sorter, 4x25 cassette capacity	9001387
Duo unit for Vega	8000043
Trays for Vega Duo, 2 x 10 cassette capacity	9001547
Hopper carousel with 6 hoppers, 480 cassette capacity	9001327
Hopper carousel with 8 hoppers, 640 cassette capacity	9001521
Vega Collection tray	9002054
Scanner Kit & Packaging	B35310050
13.3 in. All-in-one PC Package	B81300033



Epredia VEGA laser cassette printer with optional sorter system



Epredia VEGA laser cassette printer with optional Duo sorter system

For more information, contact your local Epredia Account Executive or visit epredia.com

